

Gall Duct Forceps, Probes, Scissors, Kidney Forceps, Urinary Instruments, Trocars



Gall Duct Forceps



03-5622
LOWER
180 mm



03-5623
LAHEY (SWEET)
190 mm



03-5624
LAHEY
220 mm




03-5625
DESJARDINS
Gall duct-and cystic forceps
210 mm







03-5626
NISSEN
Gall duct-and cystic forceps
215 mm

Gall Duct Forceps and Choledochus Forceps



Gall Bladder Forceps and Gall Duct Scis:



Gall Duct Probes

Ø 3.5 mm



03-5639
HAJEK
Irrigating cannula,
silver
165 mm

Ø 5 mm




03-5640
MOYNIHAN
Gall stone probe
and scoop of tin,
malleable.
365 mm

Ø 5 mm



03-5641
MAYO
Gall duct probes,
malleable
330 mm

Ø 3.5 mm



03-5642
UCHSNER
Gall duct probes, flexible
400 mm

Ø 3.8 mm



03-5643




03-5644
Ø 4.3mm



03-5645
Ø 6.8mm

03-5643 - 03-5645
DESJARDINS
Gall duct probes,
flexible
285 mm

Gall Duct Dilators & Gall Stone Scoops






03-5646 - 03-5668
DITTEL
Metal bougies and
dilating sounds, stiff
345 mm

03-5669 - 03-5681
BAKES
Common bile duct
dilators, malleable
03-5669-03-5673 315mm
03-5674-03-5681 320mm
Set of 9 pieces (**03-5671-03-5679**)
in canvas roll

03-5682 - 03-5685
DESJARDINS
Probes and catch
with flexible shaft
300 mm


03-5686 - 03-5693
LUER-KOERTE
modif with malleable shaft
310 mm

Gall Stone Forceps




Kidney Clamps


Atraumata



Kidney and Pedicle Clamps



03-5709
STILLE
220 mm



03-5710
HERRICK
225 mm



03-5711
WERTHEIM
225 mm

Kidney Elevating Forceps



03-5712
COLLIN
285 mm



03-5713
DARTIGUES
255 mm

Kidney Stone Forceps



03-5714



03-5715



03-5716



1/4
03-5717



1/2

03-5714 - 03-5717
RANDALL
225 mm



03-5718



1/4
03-5719




1/2

03-5718 - 03-5719
RANDALL
190 mm

Bladder Retractors




03-5720
JUDD-MASSON
Bladder-and suprapubic
retractor



03-5721
MILLIN
Center blade for (**05-6993**)
retro pubic prostatectomy




03-5722
MILLIN
Bladder retractor
03-5722 Spreader only
03-5723 Blades 57x25mm
03-5724 Blades 82x25mm
03-5725 Center blade 105x45mm



03-5726
LEGUEU
Bladder retractor
Center blade



03-5727
LEGUEU
Bladder spatula
260 mm




03-5728
KOCHER
Bladder spatula
235 mm



03-5729
MATHIEU
Foreign body forceps
Length of shaft
190 mm




03-5730
MILLIN
Lobe grasping
forceps
235 mm



03-5731
Ureter forceps
240 mm

Atraumata



03-5732
MILLIN
Lobe grasping
forceps
245 mm




03-5733
MILLIN
Bladder neck spreader
285 mm



03-5734
GUYON
Catheter guide
(introducer)
365 mm

Charr.
(French)


03-5735	6
03-5736	7
03-5737	8
03-5738	9
03-5739	10
03-5740	11
03-5741	12
03-5742	13
03-5743	14
03-5744	15
03-5745	16
03-5746	18
03-5747	20



03-5735 - 03-5747
Female catheters,
metal
160 mm


Charr.
(French)

03-5748	8
03-5749	9
03-5750	10
03-5751	11
03-5752	12
03-5753	13
03-5754	14
03-5755	15
03-5756	16
03-5757	17
03-5758	18
03-5759	19
03-5760	20
03-5761	21
03-5762	22
03-5763	23
03-5764	24
03-5765	25
03-5766	26
03-5767	27
03-5768	28
03-5769	29
03-5770	30



03-5748 - 03-5770
DITTEL
Urethral bougies
and dilating sounds
345 mm

03-5771	Ben. 24/Charr. 12
03-5772	Ben. 25
03-5773	Ben. 26/Charr. 13
03-5774	Ben. 27
03-5775	Ben. 28/Charr. 14
03-5776	Ben. 29
03-5777	Ben. 30/Charr. 15
03-5778	Ben. 31
03-5779	Ben. 32/Charr. 16
03-5780	Ben. 33
03-5781	Ben. 34/Charr. 17
03-5782	Ben. 35
03-5783	Ben. 36/Charr. 18
03-5784	Ben. 37
03-5785	Ben. 38/Charr. 19
03-5786	Ben. 39
03-5787	Ben. 40/Charr. 20
03-5788	Ben. 41
03-5789	Ben. 42/Charr. 21
03-5790	Ben. 43
03-5791	Ben. 44/Charr. 22
03-5792	Ben. 45
03-5793	Ben. 46/Charr. 23
03-5794	Ben. 47
03-5795	Ben. 48/Charr. 24
03-5796	Ben. 49
03-5797	Ben. 50/Charr. 25
03-5798	Ben. 51
03-5799	Ben. 52/Charr. 26
03-5800	Ben. 53
03-5801	Ben. 54/Charr. 27
03-5802	Ben. 55
03-5803	Ben. 56/Charr. 28
03-5804	Ben. 57
03-5805	Ben. 58/Charr. 29
03-5806	Ben. 59
03-5807	Ben. 60/Charr. 30




03-5771 - 03-5807
GUYON
Urethral bougies
and dilating sounds,
without filiform guides,
with thread in tip
255 mm

Genito-Urinary Instruments & Circumcision Clamps



03-5808
GUYON
Set of urethral bougies and dilating sounds, in metal case, with lift out rack and bottom rack, containing 21 urethral and dilating sounds, metal.



03-5809 - 03-5811
WINKELMANN
03-5809 ∅ 12mm
03-5810 ∅ 16mm
03-5811 ∅ 20mm

Circumcision Clamps


Assisted by two haemostats or ligatures the foreskin is pushed over the bellshaped protection - rod and the glans is placed into it. By now, the wing - nut is tightened as close as possible, so the foreskin in annular compressed. The foreskin is cut of with a scalpel, right nest to the compression ring. The circumcison clamp remains a few minutes to avoid bleeding. In case that any large vein, specially the frenular artery, is cut through, a ligature must be set. The advantage of the method described is the cosmeticy favourable scar which is left after rejection of the small necrotized portion.

Trocars

Luer Lock



03-5812
KLIMA ROSEGGER
35 X 1,8 mm
Sternal puncture needle



03-5813
Universal trocar
with 4 tips



03-5814 - 03-5818
NELSON

	Inside diam	Outside diam
03-5814	2,2 mm	3,0 mm
03-5815	3,0 mm	4,0 mm
03-5816	4,0 mm	5,0 mm
03-5817	5,0 mm	6,0 mm
03-5818	5,5 mm	6,5 mm

for hydrocele and abdominal
puncture



03-5819 - 03-5822

	Inside diam	Outside diam
03-5819	2,2 mm	3,0 mm
03-5820	3,0 mm	4,0 mm
03-5821	4,0 mm	5,0 mm
03-5822	5,0 mm	6,0 mm


for abdominal puncture,
with stop-cock and lateral
drain



03-5823 - 03-5824
BUELAU

	Inside diam	Outside diam
03-5823	6,5 mm	8,0 mm
03-5824	9,0 mm	10,0 mm

Empyema trocar



03-5825 - 03-5828
OCHSNER

	Inside diam.	Outside diam
03-5825	3,0 mm	4,0 mm
03-5826	3,6 mm	4,6 mm
03-5827	4,3 mm	5,3 mm
03-5828	5,0 mm	6,0 mm



03-5829
PIERCE
Antrum trocar
ø 4,0 mm



03-5830
COAKLEY
Antrum trocar
ø 3,0 mm



Plot # 13/B, Small Industrial Estate,
Sialkot 51310 - Pakistan.
Tel: +92 53 3552444
Fax: +92 53 3550700
Email: info@g-devices.com / glomeddevices@gmail.com
URL: www.g-devices.com