

## Plaster Instruments



## Plaster Instruments

### Plaster Knives



with ebony handle

**07-9001**  
HOPKINS  
200 mm



**07-9002**  
ESMARCH  
180 mm

### Plaster Saw



**07-9003**  
ENGEL  
160 mm

### Bandage and First-Aid scissors



1/2

**07-9004 - 07-9012** 145mm

|                |             |                |             |
|----------------|-------------|----------------|-------------|
| <b>07-9004</b> | black       | <b>07-9008</b> | yellow      |
| <b>07-9005</b> | red         | <b>07-9009</b> | white       |
| <b>07-9006</b> | blue        | <b>07-9010</b> | orange      |
| <b>07-9007</b> | light green | <b>07-9011</b> | violet      |
|                |             | <b>07-9012</b> | bright pink |




1/2

**07-9013 - 07-9021** 180mm

|                |             |                |             |
|----------------|-------------|----------------|-------------|
| <b>07-9013</b> | black       | <b>07-9017</b> | yellow      |
| <b>07-9014</b> | red         | <b>07-9018</b> | white       |
| <b>07-9015</b> | blue        | <b>07-9019</b> | orange      |
| <b>07-9016</b> | light green | <b>07-9020</b> | violet      |
|                |             | <b>07-9021</b> | bright pink |

# Plaster Scissors

For small plaster dressings






# Plaster Shears

## Special features:

The construction of the open jaws avoids tiresome clogging by plaster remnants and allows furthermore an easy cleaning. The shears do not glide off the plaster owing to the deep teething on the jaws. The knives with hollow grinding have a better cutting ability and therefore less force is required. The knives are guided by ball bearings.



**07-9029**  
STILLE  
260 mm



**07-9030**  
STILLE  
360 mm

## STILLE Plaster Shears

Blade with reinforced back



**07-9032**  
STILLE  
230 mm

**07-9033**  
Spare blade



**07-9034**  
STILLE  
370 mm

**07-9035**  
Spare blade



**07-9036**  
STILLE  
240 mm

**07-9037**  
Spare blade

## Plaster Spreader, Plaster Breakers



**07-9038**  
**HENNIG**  
270 mm



**07-9039**  
**MOD.USA**  
230 mm



**07-9040**  
290 mm



with removable jaws



**07-9041 - 07-9042**  
**WOLFF**  
**07-9041** 180 mm  
**07-9042** 240 mm  
plaster breakers



Plot # 13/B, Small Industrial Estate,  
Sialkot 51310 - Pakistan.  
Tel: +92 53 3552444  
Fax: +92 53 3550700  
Email: [info@g-devices.com](mailto:info@g-devices.com) / [glomeddevices@gmail.com](mailto:glomeddevices@gmail.com)  
URL: [www.g-devices.com](http://www.g-devices.com)