


Suction Instruments



Vacuum Curettes for Aspiration Lipectomy Suction Instruments



07-9043 - 07-9048

BIEMER

- 07-9043 ∅ 4mm
- 07-9044 ∅ 6mm
- 07-9045 ∅ 8mm
- 07-9046 ∅ 10mm
- 07-9047 ∅ 12mm
- 07-9048 ∅ 14mm

Curettes for aspiration lipectomy, graduated, slightly curved

07-9049 - 07-9050

BIEMER

- 07-9049 ∅ 4 mm
 - 07-9050 ∅ 6 mm
- Curettes for aspiration lipectomy, with three apertures, slightly curved malleable

07-9051 - 07-9053

BIEMER

- 07-9051 ∅ 8 mm
 - 07-9052 ∅ 10 mm
 - 07-9053 ∅ 12 mm
- Curettes for aspiration lipectomy, with three apertures, slightly curved malleable

07-9054

For aspirating fluid from the excretory ducts of large glands and for keeping the operating area clear 170 mm

07-9055

Trocar
For puncturing 150 mm

07-9056

For clearing large quantities of fluid 190 mm

Suction Instruments

07-9057 150mm

Suction handle, with vacuum regulator and finger cut-off, brass, mirror-finish chromium-plated



olive
Ø 9 mm

07-9058 170mm

Suction handle, with hand regulation valve and quick-shutter, brass, mirror-finish chromium-plated



★ Both suction handles have bayonet fitting olive and can be used for the suction cannulae. Ø 9 mm



07-9059
Trocár
For puncturing
170 mm



07-9060
Hollow needle,
sharp
150 mm



07-9061
double curved
For general
surgery
230 mm




07-9062
For tonsillectomy,
oral and laryngeal
operations
190 mm



07-9063
For tonsillectomy,
oral and laryngeal
operations
220 mm

★ With bayonet fitting to be used with suction handles (07-9057 - 07-9058)

Suction Instruments



07-9064
For tonsillectomy,
oral-and laryngeal
operations
180 mm

07-9065
Suction cannula
for tonsillectomy,
oral-and laryngeal
operations
220 mm


07-9066
Suction cannula
for tonsillectomy,
oral-and laryngeal
operations
180 mm

07-9067
S-shaped blunt cannula
for aspirating fluid from
the excretory ducts of
large glands and for
keeping the operating
area clear
170 mm

07-9068
POOL
Suction cannula,
with outer tube
to be screwed of
225 mm

07-9069
POOL
Suction cannula,
with outer tube
to be screwed off
225 mm

Suction Instruments




Suction Instruments



07-9076 Ø 1,5 mm
07-9077 Ø 2,0 mm
07-9078 Ø 2,5 mm
07-9079 Ø 3,0 mm

07-9076 - 07-9079
PLESTER
with finger cut-off and
stylet, thumb piece
parallel to curve, with
LUER hub
195 mm




07-9080 - 07-9084
FERGUSON
with finger cut-off and
stylet, thumb piece
90° to curve, with
LUER hub
180 mm



07-9085 Ø 4mm
FERGUSON
insulated, with finger cut-off
and stylet, thumb piece 90°
to curve, with LUER hub
180 mm



07-9086 Ø 3mm
07-9087 Ø 5mm
FRAZIER
with finger
cut-off and stylet, thumb
piece 90° to curve,
with LUER hub
180 mm



07-9088
5,0 mm diam., Tip
6,5 mm wide, Nasal
suction tube, with
LUER hub
180 mm



07-9089 4,5 mm
07-9090 5,5 mm
Nasal suction
tubes
175 mm



07-9091
GORNEY
Suction dissector
180 mm



07-9092 4 mm wide
07-9093 5 mm wide
Suction raspatories,
sharp, bayonet-shaped,
with stylet
200 mm



07-9094
LOPEZ-REINKE
Tonsil dissector with
suction device
235 mm

Suction Instruments



"Supporting Healthcare Globally"

07-9095
WULLSTEIN
Suction tube, conical, silver, olive
with LUER hub
160 mm



07-9096 1,5 mm
07-9097 2,0 mm
WULLSTEIN
Suction tube, cylindrical, silver, with
stylet olive with LUER hub
160 mm



Usable length

YASARGIL usable length diam.

07-9098	130 mm	1,5 mm
07-9099	150 mm	1,5 mm
07-9100	180 mm	1,5 mm
07-9101	220 mm	1,5 mm

YASARGIL usable length diam.

07-9106	130 mm	2,5 mm
07-9107	150 mm	2,5 mm
07-9108	180 mm	2,5 mm
07-9109	220 mm	2,5 mm


07-9102	130 mm	2,0 mm
07-9103	150 mm	2,0 mm
07-9104	180 mm	2,0 mm
07-9105	220 mm	2,0 mm

07-9110	130 mm	3,5 mm
07-9111	150 mm	3,5 mm
07-9112	180 mm	3,5 mm
07-9113	220 mm	3,5 mm
07-9114	180 mm	4,5 mm

olive
ø 6,5-9 mm
with LUER hub

07-9115 - 07-9117
LANDOLT
150 MM
diam ø 3 mm

07-9115 8 mm / curved
07-9116 11mm / curved
07-9117 15mm / curved







07-9118 500 mm
07-9119 350 mm
with rubber tip,
with LUER hub

07-9120 2,0 mm diam
07-9121 2,5 mm diam
07-9122 3,0 mm diam
with finger cut-off,
thumb piece 90° to
curve, with LUER hub
250 mm



Combined Suction and Irrigation Cannula Coagulation Suction Cannulae





Plot # 13/B, Small Industrial Estate,
Sialkot 51310 - Pakistan.
Tel: +92 53 3552444
Fax: +92 53 3550700
Email: info@g-devices.com / glomeddevices@gmail.com
URL: www.g-devices.com