


[illegible]

Wire Extension



06-8554 - 06-8560
KIRSCHNER
Wire traction bows,
to be used with
KIRSCHNER-wires
up to 2.5 mm diam.,
With 3 hooks,
stainless steel

06-8554
06-8555
06-8556
06-8557
06-8558
06-8559
06-8560

Inside dimensions
width

200 mm
155 mm
135 mm
120 mm
105 mm
95 mm
95 mm

depth

155 mm
155 mm
130 mm
120 mm
100 mm
110 mm
70 mm

for the hip
for the femur
for the knee
for the knee
for the knee
for the calcaneus
for the arm and elbow



06-8561
Wrench, double - ended,
for wire traction bows



06-8562
Wire tightener, with
clamping jaws



06-8563
Traction hook only,
large size




06-8564
Traction hook only,
small size




06-8565
Fixation clamp, to secure clamping
jaws of **06-8562** stainless steel


Wire Extension




06-8566
Wire traction bow, adjustable,
with 3 hooks, small size



06-8567
Wire traction bow, adjustable,
with 3 hooks, large size



06-8568
Wire traction bow, adjustable,
with 3 hooks, small size for
elbow and calcaneus, for
adolescents



06-8569 max. 105 mm
06-8570 max. 160 mm
06-8571 max. 200 mm

06-8569 80-95 mm
06-8570 90-120 mm
06-8571 140-170 mm

dial upto
06-8569 1,6 mm
06-8570 2,0 mm
06-8571 2,5 mm

American Type

06-8572
Fixation disc, slotted

06-8573
idem, aluminum, transparent
to X-rays

06-8569 - 06-8571
KIRSCHNER
Wire traction bows, modified,
adjustable in width, holds wire
securely, with 1 hook

06-8574
Fixation disc, slotted

06-8575
Screw driver, for above
fixations discs

Traction Bows, Distractor




Boehler traction bows, to be used with STEINMANN pins, 4mm diam

Boehler traction bows, to be used with STEINMANN pins, 5mm diam

		Inside width	For
06-8576 with pin (09-8588)		16x9cm	Arm
06-8577 without pin	06-8582 without pin		
06-8578 with pin (09-8590)		21x11cm	Femur and knee
06-8579 without pin	06-8583 without pin		
06-8580 with pin (09-8592)		21x15cm	Femur and knee
06-8581 without pin			



06-8599
SCHWARZ
Traction bow for fingers, to be used with KIRSCHNER-wire up to 1.6 mm diam.



06-8600
Small distractor as a repositioning device for distraction and compression of small bones, for wires with 3.0 mm diam.



06-8584 - 06-8598
STEINMANN
triangular shaft ended and trocar point

diam ϕ	4,0 mm	5,0 mm
Length		
80 mm	06-8584	
100 mm	06-8585	
120 mm	06-8586	
140 mm	06-8587	
150 mm	06-8588	06-8594
160 mm	06-8589	
180 mm	06-8590	06-8595
200 mm	06-8591	06-8596
210 mm	06-8592	
250 mm	06-8593	06-8597
300 mm		06-8598



06-8601
Protection cap for 4,0 mm STEINMANN Pins



06-8602
Protection cap for 5,0 mm STEINMANN Pins



1/2

06-8603
REILL
Drilling chuck,
cannulated,
for shanks from
2.0 - 2.5 mm diam
115 mm




1/2

06-8604
Drilling chuck,
cannulated,
for shanks from
0.8 - 2.0 mm diam
150 mm



1/2

06-8605
Drilling chuck,
cannulated,
for shanks up to
4.5 mm diam
with key
140 mm



1/2

06-8606
LOUTE
Wire tightener, improved pattern,
with movable tip,
for wires up to 1.0 mm diam.
220 mm

Chuck jaws made of material
used for implants



1/2

06-8607
Drilling chuck, cannulated,
for shanks up to 6.35 mm diam.,
To be tightened by hand
100 mm



1/2

06-8608
Drilling chuck, cannulated,
for shanks up to 4.5 mm diam.,
And for triangular shaft width
across flats 4.7 mm (schanz screws,
steinmann pins). with key **06-8609**

06-8609




1/2


06-8610
Drilling chuck, cannulated,
for shanks up to 6.35 mm diam.,
With key **(06-8611)**
160 mm

06-8611

Cerclage Instruments



06-8612
275 mm




06-8613
285 mm



06-8614
305 mm




06-8612 - 06-8614
DEMEL
Wire guides



06-8615
240 mm




06-8616
270 mm



06-8617
225 mm

JORDAN instruments for subcutaneous Tibia Cerclage



06-8618
JORDAN
Wire tightener



06-8619
JORDAN
Wire twister
170 mm



06-8620
JORDAN
Wire catcher
180 mm



06-8621
JORDAN
Probe, cannulated
170 mm

Wire Holding forceps, Wire Tightening Pliers Flat-nosed Pliers DUROGRIP with Hard N Inserts




06-8622
DUROGRIP
Holding forceps for cerclage
wire, with lateral jaws to avoid
wire slipping off
170 mm



06-8623
DUROGRIP
Wire, tightening pliers,
Ligature pliers
150 mm



06-8624
DUROGRIP
Wire, tightening pliers,
Ligature pliers
150 mm



06-8625
with serrated jaws
140 mm



06-8626
one jaw smooth,
one jaw serrated
140 mm




06-8627
with serrated jaws
06-8628
with longitudinal and
transversal grooves
170 mm



06-8629
190 mm

Wire Holding forceps, Flat-nosed Pliers, Pliers for bending and cutting wires

GLOMEDdevices
"Supporting Healthcare Globally"




06-8630
130 mm

06-8631
180 mm

06-8630 - 06-8631
Extraction pliers
for boring wires



06-8632
Flat-nosed pliers,
parallel
185 mm




06-8633
Vise grip
200 mm




06-8634
Round Pliers
180mm

06-8634 - 06-8637



06-8635 145 mm
06-8636 180 mm
06-8637 180 mm
Extraction pliers for wires



06-8638
Bending pliers
for wires up to
1.8 mm diam.


Wire Cutters and wire Cutting Scissors DUROTIP with Carbide Cutting Edges



06-8639 **06-8640**
ANGLE DUROTIP
130 mm ANGLE
130 mm

for soft wire
up to diam.
for hard wire
up to diam.


0,5 mm 0,6 mm
0,3 mm 0,4 mm



06-8641 **06-8642**
120 mm DUROTIP
120 mm

for soft wire
up to diam.
for hard wire
up to diam.

1,2 mm 1,2 mm
0,8 mm 0,8 mm



06-8643
HARVEY
125 mm

0,4 mm
0,3 mm



06-8644
150 mm

for soft wire
up to diam
for hard wire
up to diam

1,5 mm
1,0 mm




06-8645
160 mm

1,5 mm
1,2 mm



06-8646
140 mm
DUROTIP
1,5 mm

1,0 mm



06-8647
180 mm

for soft wire
up to diam. 2,0 mm
For hard wire 1,5 mm
up to diam
specially hardened
cutting edges, front
and lateral cutting
action



06-8648
REILL
175 mm

for soft wire
up to diam. 2,0 mm
For hard wire 1,6 mm
up to diam
specially hardened
cutting edges, for
cutting off wires close
to the bone



DUROTIP

06-8649
235 mm

for soft wire
up to diam. 3,0 mm
For hard wire 2,2 mm
up to diam

Wire Cutting Pliers, Lateral and Front Cutting Action



06-8650
240 mm

for soft wire
up to diam. 2,0 mm
For hard wire
up to diam 1,6 mm



06-8651
145 mm

for soft wire
up to diam. 1,5 mm
For hard wire
up to diam 1,0 mm



06-8652
145 mm

for soft wire
up to diam. 2,0 mm
For hard wire
up to diam 1,5 mm



06-8653
220 mm

for soft wire
up to diam. 3,5 mm
For hard wire
up to diam 2,5 mm

Wire Cutting Pliers



06-8654
320 mm

for soft wire
up to diam. 3,0 mm
For hard wire
up to diam. 3,0 mm




Cutter for cutting off
KIRSCHNER wires,
Tibia head screws,
bolts, pins, up to 3 mm diam.



DUROTIP

06-8655
Cutter for cutting off
KIRSCHNER wires,
Tibia head screws,
bolts, pins, up to 6 mm diam.
470 mm


Finger Nail Instruments



06-8656
Nail extracting
forceps
135 mm




06-8657
Nail extracting
forceps
145 mm



06-8658
Finger nail drill
with 3 drill points

Drill points only



06-8662
Nail splitting
forceps
135 mm



06-8663
SYSTRUNK
Scissors, for splitting
finger nails
130 mm



06-8664
Finger ring saw
150 mm

Spare blade
06-8665



Plot # 13/B, Small Industrial Estate,
Sialkot 51310 - Pakistan.
Tel: +92 53 3552444
Fax: +92 53 3550700
Email: info@g-devices.com / glomeddevices@gmail.com
URL: www.g-devices.com