



GLOMEDevices

"Supporting Healthcare Globally"




Dental
Instruments
Catalogue






Impression Trays

Made of 18/8 Stainless Steel, Solid, Set of 10




All The Above Impression Trays are Available in These 5 Sizes



UPPER		U1	U2	U3	U4	U5
	A	62	60	53	54	52
	B	83	75	73	69	68


All Dimension are in mm



LOWER		L1	L2	L3	L4	L5
	A	63	61	57	60	53
	B	82	81	78	74	71

All Dimension are in mm


Perforated, Set of 10






Impression Trays

Made of 18/8 Stainless Steel, Solid Full Denture Set of 8




All The Above Impression Trays are Available in These 4 Sizes




UPPER		XL	L	M	S
	A	70	67	64	60
	B	79	76	70	73

All Dimension are in mm



LOWER		XL	L	M	S
	A	64	64	60	57
	B	86	70	73	70

All Dimension are in mm




Perforated Full Denture Set of 8




Impression Trays

Made of 18/8 Stainless Steel, Solid Set of 6



All The Above Impression Trays are Available in These 4 Sizes





UPPER		U1	U2	U3
	A	67	64	57
	B	81	71	70

All Dimension are in mm

LOWER		L1	L2	L3
	A	60	57	54
	B	81	79	76

All Dimension are in mm




Perforated Set of 6



Impression Trays


Made of 18/8 Stainless Steel, English Pattern Set of 8




2889
Lower



2890
Upper



2891
Lower



2892
Upper




2893
Lower



2894
Upper

All The Above Impression Trays are Available in These 4 Sizes




UPPER	Size	A	B	C
	0	67	50	12
	1	72	51	13
	2	78	55	15
	3	80	58	15

All Dimension are in mm

LOWER	Size	A	B	C
	0	70	50	14
	1	77	51	15
	2	79	55	15
	3	80	60	17

All Dimension are in mm






Impression Trays

Made of 18/8 Stainless Steel, Solid Full Denture Set of 8




All The Above Impression Trays are Available in These 4 Sizes




UPPER		XL	L	M	S
	A	70	67	64	60
	B	79	76	70	73

All Dimension are in mm



LOWER		XL	L	M	S
	A	64	64	60	57
	B	86	70	73	70

All Dimension are in mm



Made of 18/8 Stainless Steel, Perforated Full Denture Set of 8



Impression Trays

Made of 18/8 Stainless Steel, Solid Full Denture Set of 8



2911
U0



2912
U1



2913
U2



2914
U3



2915
L0



2916
L1



2917
L2




2918
L3

Upper Sizes: 0,1,2,3


Depth of Tray: 15,9 mm, Height of Plate: 12,5 mm

Lower Sizes: 0,1,2,3


Depth of Tray: 12,9 mm




2919
U0



2920
U1



2921
U2




2922
U3




2923
L0



2924
L1



2925
L2



2926
L3

Made of 18/8 Stainless Steel, Perforated Full Denture Set of 8



Impression Trays

Impression Trays For Crown & Bridge Work



2927
U1P



2928
L1P

SOLID



2929
U1P



2930
L1P

PERFORATED



2932
Adjustable
Perforated



2933
Adjustable
Solid

PERFORATED

SOLID



2931-A
For upper right jaw
For lower left jaw

2931-B
For upper left jaw
For lower right jaw

Available in 4 sizes

- 1.....44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm




2931-C
For upper right jaw
For lower left jaw

2931-D
For upper left jaw
For lower right jaw

Available in 4 sizes

- 1.....44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm



Factory Address:

Plot # 13/B, Small Industrial Estate,
Opposite to Halifax & Co (Pvt) Ltd.
Sialkot 51310 - Pakistan.

Tel: +92 52 3552444 / Fax: +92 52 3550700

Email: info@g-devices.com / glomedevices@gmail.com

URL: www.g-devices.com



ISO 9001:2008



ISO 13485:2003

