




Disposable
Instruments

www.g-devices.com

Usage Instructions

Glomed provides High Quality Medical instruments for the whole world market that are backed by our 35 years of experience in this field. The quality of our Medical Instruments is supreme and we invite you to use them with the confidence and knowledge that you are using the finest premium products available.

These instruments are delivered sterile and ready to use.

 Instruments must be used for their specified purpose and incorrect usage could damage the instrument. Our instruments are designed to be used by competent surgeons, who have good knowledge of their feature and how they should be used. Any other use can compromise the safety of the user and the patient. It is the responsibility of the surgeon to choose the most suitable instrument for the surgical technique being performed, based on his experience and expertise.

 Only to be used by personnel with the required specialist knowledge and training for the purpose intended.

STERILIZATION:

100% ETO CYCLES

Concentration ETO: 850 ± 50mg/l

Temperature: 37°C - 47°C

Exposure Time: 3-4 hours

Humidity: 70% RH minimum

Along with these steps mentioned Gamma irradiation protocol can also be used with these products. Above mentioned sterilization cycles represent industry standard and should be capable of producing a sterile device. Due to variations in sterilization equipment and device bioburden in clinical use,, Eye technology is not able to produce specific cycle parameters. It is the responsibility of each user to perform a validation and verification of the sterilization cycle to ensure an adequate sterility assurance level of our products.

SHELF LIFE

3 years for Sterile Packaging

RECYCLING:

Glomed recommends recycling of the instruments after use. When reprocessing medical instruments please follow your local Government Health and Safety procedures.



TABLE OF CONTENTS



Scissors
01 - 07



Retractors
18 - 28



Forceps
08 - 09



Suture
29 - 30



Artery Forceps
10 - 15



Bone Surgery
31




Cotton Swab Forceps
16 - 17

GLOMED DEVICES


SCISSORS

EN Operating Scissors
DE Chirurgische Scheren
ES Tijeras Quirúrgicas



Standard

04-0010-11	11.5 cm	04-0011-11
04-0010-13	13 cm	04-0011-13
04-0010-14	14.5 cm	04-0011-14
04-0010-16	16.5 cm	04-0011-16
04-0010-18	18.5 cm	04-0011-18
04-0010-20	20 cm	04-0011-20



Standard

04-0012-11	11.5 cm	04-0013-11
04-0012-13	13 cm	04-0013-13
04-0012-14	14.5 cm	04-0013-14
04-0012-16	16.5 cm	04-0013-16
04-0012-18	18.5 cm	04-0013-18
04-0012-20	20 cm	04-0013-20



Standard

04-0014-11	11.5 cm	04-0015-11
04-0014-13	13 cm	04-0015-13
04-0014-14	14.5 cm	04-0015-14
04-0014-16	16.5 cm	04-0015-16
04-0014-18	18.5 cm	04-0015-18
04-0014-20	20 cm	04-0015-20

Operating Scissors **EN**
 Chirurgische Scheren **DE**
 Tijeras Quirúrgicas **ES**

GLOMED DEVICES

SCISSORS



Incision
 04-0040-14
 14 cm




Deaver
 04-0050-14 14 cm 04-0051-14




Mayo
 04-0060-14 14.5 cm 04-0061-14
 04-0060-15 15 cm 04-0061-15
 04-0060-17 17 cm 04-0061-17
 04-0060-23 23 cm 04-0061-23



Mayo-Harrington
 04-0070-22 22.5 cm 04-0071-22
 04-0070-30 30 cm 04-0071-30



Mayo-Stille
 04-0080-15 15 cm 04-0081-15
 04-0080-17 17 cm 04-0081-17




Mayo-Stille
 04-0090-16 16 cm 04-0091-16
 04-0090-19 19 cm 04-0091-19
 04-0090-21 21 cm 04-0091-21

GLOMED DEVICES

SCISSORS

EN Operating Scissors, Dissecting Scissors
 DE Chirurgische Scheren, Präparierscheren
 ES Tijeras Quirúrgicas, Tijeras De Quirúrgicas




Metzenbaum

04-0160-11 11.5 cm 04-0161-11



Metzenbaum

04-0170-15 15.5 cm 04-0171-15



Kelly

04-0180-16 16 cm 04-0181-16



Metzenbaum
14.5 cm



04-0190-14



04-0192-14



04-0194-14



04-0191-14



04-0193-14



04-0195-14



Dissecting Scissors EN

Präparierscheren DE

Tijeras Quirúrgicas ES

GLOMED DEVICES

SCISSORS

**Metzenbaum-Nelson**


04-0200-18	18 cm	04-0201-18
04-0200-20	20 cm	04-0201-20
04-0200-23	23 cm	04-0201-23
04-0200-25	25 cm	04-0201-25
04-0200-28	28 cm	04-0201-28
04-0200-30	30 cm	04-0201-30

**Metzenbaum-Nelson**

04-0210-18	04-0210-20
18 cm	20.5 cm

**Metzenbaum-Nelson**

04-0220-14	14.5 cm	04-0221-14
------------	---------	------------

**Lexer**

04-0230-16	16 cm	04-0231-16
------------	-------	------------

**Fomon**

04-0241-14
14 cm

**Lexer-Baby**

04-0250-10	10 cm	04-0251-10
------------	-------	------------

GLOMED DEVICES

SCISSORS

EN Episiotomy Scissors & Umbilical Scissors

DE Episiotomiescheren & Nabelschnurscheren

ES Tijeras Para Episitomia & Tijeras Umbilicales



Modell USA
04-0450-10
10.5 cm




Busch
04-0470-16
16cm

Operating and Gum Scissors **EN**
 Chirurgische und Zahnfleisch Scheren **DE**
 Tijeras Quirurgicas Para encías **ES**

GLOMED DEVICES

SCISSORS



Wagner

04-0490-12 12 cm 04-0491-12



La Grange

04-0500-11
11.5 cm




Iris

04-0520-09 9 cm 04-0521-09
 04-0520-10 10 cm 04-0521-10
 04-0520-11 11.5 cm 04-0521-11



Stevens
11 cm



04-0530-11

04-0531-11

04-0540-11

04-0541-11

GLOMED DEVICES

SCISSORS

EN Stitch Scissors, Wire Cutting Scissors

DE Fadenscheren, Drahtschere

ES Tijeras para hilos, Tijeras para alambre




Spencer
04-0600-13
13 cm



Littauer
04-0610-14
14 cm



Beebee



04-0630-10 10 cm 04-0631-10



04-0640-12 12 cm 04-0641-12

GLOMED DEVICES

FORCEPS

EN Dressing Forceps
DE Anatomische Pinzetten
ES Pinzas de disección



10.5 cm 05-0010-10
11.5 cm 05-0010-11
13 cm 05-0010-13
14.5 cm 05-0010-14
16 cm 05-0010-16
18 cm 05-0010-18
20 cm 05-0010-20
23 cm 05-0010-23
25 cm 05-0010-25
30 cm 05-0010-30



11.5 cm 05-0020-11
13 cm 05-0020-13
14.5 cm 05-0020-14
16 cm 05-0020-16
18 cm 05-0020-18
20 cm 05-0020-20
25 cm 05-0020-25
30 cm 05-0020-30



05-0040-12
12 cm



05-0040-15
15 cm

Dissecting Forceps, Jeweller Type Forceps **EN**
 Chirurgische Pinzetten, Uhrmacher-Pinzetten **DE**
 Pinzas de disección con dientes, Pinzas de relojero **ES**

GLOMED DEVICES FORCEPS



	1:2
10.5 cm	05-0400-10
11.5 cm	05-0400-11
13 cm	05-0400-13
14.5 cm	05-0400-14
16 cm	05-0400-16
18 cm	05-0400-18
20 cm	05-0400-20
23 cm	05-0400-23
25 cm	05-0400-25
30 cm	05-0400-30



	2:3
11.5 cm	05-0410-11
13 cm	05-0410-13
14.5 cm	05-0410-14
16 cm	05-0410-16
18 cm	05-0410-18
20 cm	05-0410-20
25 cm	05-0410-25



Adson

05-0500-12
12 cm05-0500-15
15 cm

Micro-Adson

05-0630-12
12 cm05-0630-15
15 cm

Micro-Adson

05-0640-12
12 cm05-0640-15
15 cm

Feilchenfeld

05-0780-07
7.5 cm
05-0780-09
9 cm
05-0780-11
11.5 cm



GLOMED DEVICES

ARTERY FORCEPS

EN Haemostatic Forceps, Organ Grasping Forceps

DE Arterienklemmen, Organfasszangen

ES Pinzas hemostáticas, Pinzas para sacar órganos




Dissecting and Ligature Forceps, Haemostatic Forceps EN

Präparierklemmen, Arterienklemmen DE


Pinzas para preparación y ligadura, Pinzas hemostáticas ES

GLOMED DEVICES

ARTERY FORCEPS

**Mixter-Baby**06-0150-01
14 cm06-0150-02
13.5 cm**Hartmann**
10 cm

06-0160-10 06-0161-10

**Hartmann**
10 cm

06-0170-10 06-0171-10

**Micro-Mosquito**

06-0180-10 10 cm 06-0181-10
06-0180-12 12.5 cm 06-0181-12
06-0190-10 10 cm 06-0191-10
06-0190-12 12.5 cm 06-0191-12

**Halsted-Mosquito**

06-0200-12 12.5 cm 06-0201-12
06-0200-14 14 cm 06-0201-14
06-0210-12 12.5 cm 06-0211-12
06-0210-14 14 cm 06-0211-14

GLOMED DEVICES

ARTERY FORCEPS

EN Haemostatic Forceps

DE Arterienklemmen

ES Pinzas Hemostáticas



06-0240-14 06-0241-14

Providence-Hospital
14 cm



06-0270-14 06-0271-14

Kelly
14 cm




06-0280-16 06-0281-16

Kelly-Rankin
16 cm



06-0300-14 06-0301-14 06-0310-14 06-0311-14

Crile
14 cm



06-0320-16 06-0321-16

Crile-Rankin
16 cm


Dissecting and Ligature Forceps, Haemostatic Forceps EN

Präparierklemmen, Arterienklemmen DE

Pinzas para preparación y ligadura, Pinzas hemostáticas ES

GLOMED DEVICES

ARTERY FORCEPS




Adson
18 cm

06-0360-18



06-0361-18



Rochester-Carmalt

06-0370-16

16 cm

06-0370-20

20 cm



06-0371-16

06-0371-20



Kocher

Delicate06-0390-16
16 cm

GLOMED DEVICES

ARTERY FORCEPS

EN Dissecting and Ligature Forceps, Haemostatic Forceps

DE Präparierklemmen, Arterienklemmen

ES Pinzas para preparación y ligadura, Pinzas hemostáticas



Bengolea



06-0450-20	20 cm	06-0451-20
06-0450-22	22 cm	06-0451-22
06-0450-24	24 cm	06-0451-24
06-0450-26	26 cm	06-0451-26



06-0460-20	20 cm	06-0461-20
06-0460-22	22 cm	06-0461-22
06-0460-24	24 cm	06-0461-24
06-0460-26	26 cm	06-0461-26



Mixter -Baby



06-0471-14	14 cm	06-0471-18
		18 cm



Adson-Baby



06-0481-14	14 cm	06-0481-18
		18 cm

Towel Clamps , Tubing Clamps **EN**
 Tuchklemmen , Schlauchklemmen **DE**
 Pinzas de campo , Pinzas para Tubos **ES**

GLOMED DEVICES

ARTERY FORCEPS



06-1100-14
14 cm



06-1110-16
16.5 cm

06-1110-18
18 cm



06-1120-15
15 cm

06-1120-18
18 cm

06-1120-20
20 cm




GLOMED DEVICES

COTTON SWAB FORCEPS

EN Dressing Forceps
DE Komzangen
ES Pinzas para vendajes



Lister
07-0030-13
13 cm



07-0070-26



Bozemann
26 cm




07-0071-26

Sponge Forceps **EN**
 Schwammzangen **DE**
 Pinzas para tapones **ES**

GLOMED DEVICES

COTTON SWAB FORCEPS



Foerster




07-0090-18 18 cm
 07-0090-20 20 cm
 07-0090-25 25 cm



07-0100-18
 07-0100-20
 07-0100-25



07-0091-18 18 cm
 07-0091-20 20 cm
 07-0091-25 25 cm



07-0111-18
 07-0111-20
 07-0111-25



Rampley



07-0130-18
 18 cm



07-0130-25
 25 cm

GLOMED DEVICES

RETRACTOR

EN Hooklets
DE Wundhäkchen
ES Ganchitos



Knapp
16.5 cm



08-0050-16
Sharp



08-0051-16
Blunt



16.5 cm



08-0060-01



sharp



08-0060-02



08-0060-03



08-0060-04



16.5 cm



08-0061-01



blunt



08-0061-02



08-0061-03



08-0061-04



16.5 cm



08-0070-01



08-0070-02



sharp



08-0070-03



08-0070-04



16.5 cm



08-0071-01



08-0071-02



blunt



08-0071-03




08-0071-04


Hooklets, Saddle Hooks **EN**
 Wundhäkchen, Sattelhäkchen **DE**
 Ganchitos, Ganchos en forma de silla **ES**

GLOMED DEVICES

RETRACTOR



Desmarres

Langenbeck-Mini
16 cmMeyerding
18 cmSenn-Green
15 cm

GLOMED DEVICES

RETRACTOR

EN Saddle Hooks, Hooklets
DE Sattelhäkchen, Wundhäkchen
ES Ganchos en forma de silla, Ganchitos



Cushing



20 cm
24 cm

8 mm



08-0120-08
08-0130-08

10 mm



08-0120-10
08-0130-10

12 mm



08-0120-12
08-0130-12

14 mm



20 cm
24 cm

08-0120-14
08-0130-14

16 mm



08-0120-16
08-0130-16

18 mm



08-0120-18
08-0130-18



Ragnell-Davis
08-0150-14
14 cm



8 x 4 mm

15 x 5.5 mm



Retractors EN

Wundhaken DE

Separadores ES

GLOMED DEVICES

RETRACTOR



Volkmann
21.5 cm



Sharp



Blunt

08-0180-01




08-0180-02



08-0180-03



08-0180-04



08-0180-06



08-0181-01



08-0181-02



08-0181-03



08-0181-04



08-0181-06



Kocher
22 cm



08-0190-01
Sharp



08-0190-02




08-0190-03



08-0190-04



08-0190-06



08-0190-08



08-0191-01
Blunt



08-0191-02




08-0191-03



08-0191-04



08-0191-06




08-0191-08

GLOMED DEVICES

RETRACTOR

EN Retractors
DE Wundhaken
ES Separadores



Volkman USA Model
22.5 cm



Volkman
21.5cm



08-0220-21
sharp



08-0221-21
blunt



Volkman
20.0cm



08-0230-20
sharp



08-0231-20
blunt




08-0240-24
24.5 cm


Retractors **EN**
 Wundhaken **DE**
 Separadores **ES**

GLOMED DEVICES

RETRACTOR




Langenbeck
 08-0310-22
 22.5 cm



Langenbeck
 08-0320-21
 21.5 cm




Lahey
 08-0330-19
 19 cm



Kocher
 23 cm

08-0350-18
 40 x 18 mm
 08-0350-20
 60 x 20 mm
 08-0350-25
 60 x 25 mm
 08-0350-30
 75 x 30 mm
 08-0350-40
 75 x 40 mm




Hosel
 26 cm

08-0440-01
 60 x 20 mm
 08-0440-02
 80 x 30 mm
 08-0440-03
 100 x 30 mm
 08-0440-04
 120 x 30 mm
 08-0440-05
 140 x 30 mm
 08-0440-06
 140 x 40 mm

GLOMED DEVICES

RETRACTOR

EN Retractors
DE Wundhaken
ES Separadores



Retractors EN

Wundhaken DE

Separadores ES

GLOMED DEVICES

RETRACTOR

08-0640-01
95 x 35 mm

08-0640-02
120 x 50 mm

08-0640-03
150 x 50 mm



Mikulicz
26 cm

08-0650-26
90 x 50 mm



Mikulicz
26 cm

08-0660-01
190 x 38 mm

08-0660-02
190 x 57 mm



Kelly
27 cm

30 x 10 mm
26 x 10 mm



Farabeuf
08-0710-12
12 cm

27 x 16 mm
23 x 16 mm



Farabeuf
08-0710-15
15 cm

30 x 13 mm
34 x 13 mm

26 x 16 mm
30 x 16 mm



Roux-Mini
08-0720-14
14 cm



Roux-Baby
08-0730-13
13 cm

22 x 15 mm
19 x 15 mm



Parker
08-0740-01
13.5 cm

GLOMED DEVICES

RETRACTOR

EN Wound Spreaders, Self-retaining

DE Wundspreizer

ES Separadores autoestáticos



Gelpi



08-1140-09

9 cm

08-1140-11

11 cm

08-1140-13

13 cm



Gelpi

08-1140-16



GLOMED DEVICES

SUTURE

EN Needle Holders

DE Nadelhalter

ES Pora-agujas




Needle Holders **EN**Nadelhalter **DE**Porta-agujas **ES**

GLOMED DEVICES

SUTURE

**Mayo-Hegar**10-0200-13
13 cm10-0200-14
14 cm10-0200-16
16 cm10-0200-18
18 cm10-0200-20
20 cm10-0200-24
24 cm10-0200-26
26 cm10-0200-30
30 cm

Delicate

**Mayo-Hegar**10-0210-16
16 cm
10-0210-18
18 cm
10-0210-20
20 cm**Olsen-Hegar**10-0290-14
14 cm
10-0290-16
16.5 cm**Gillies**
10-0300-16
16 cm




GLOMED DEVICES

BONE SURGERY

EN Bone Currettes, Double-ended

DE Scharfe Doppellöffel

ES Cucharillas de doble extremo





 Plot # 13/B, Small Industrial Estate,
Sialkot 51340 - Pakistan.

 +92. 52. 3552444

 +92. 52. 3550700

 info@g-devices.com
glomedevices@gmail.com

 www.g-devices.com

